



Technical Data Sheet

Makrite N95 Particulate Respirator-12x20 Item # MAK-9500N95

SAP Material Code & Color

| SAP Material Code | Prod REF Code | Color |
|-------------------|---------------|-------|
| MAK-9500N95 | 9500-N95 | White |

Country of Origin

China

Intended Use

Occupational (including medical) use. Reduces wearer's exposure to particles including small particle aerosols and large droplets (all non-oil aerosols). Protects from exposure to airborne particles. In a healthcare setting, protects from exposure to biological aerosols including viruses and bacteria.

Material

- Body composition
Outer material: Cotton
Filter media: Melt-blown Polypropylene
Inner material: Cotton
- Elastic Band composition
Polyester and Spandex
- Nose-piece composition
Aluminum

Dimensions

- Body Size
Length: 144 mm; Width: 117 mm; Height: 60 mm
- Length of Nose-piece
89 mm
- Length of Elastic Band
Upper: 305mm; Lower: 279mm

NIOSH Test

NIOSH approved for at least 95 percent filtration efficiency against certain non-oil based particles

Sterility

Non-sterile

Usage

Disposable, Single-use only

Packaging

20 Particulate Respirators per box, 12 boxes per case
(Total: 240 masks per case)

Recommended Storage Conditions

Store under cool, clean, and dry conditions. Avoid excessive heat
(over 40°C or 104°F)

Shelf Life

3 years

This technical data sheet is issued by A.R. Medicom Inc (Asia) Ltd