



Technical Data Sheet

Black Earloop Mask Items# VITXMK2792

SAP Material Code & Color

SAP Material Code	Prod Ref Code	Color
VITCMK2792-0100050	VITCMK2792	Black
VITGMK2792-0100050	VITGMK2792	Black

Intended Use

The procedure masks are intended to be worn by personnel to reduce inhalation of fine particulate matter including dust, pollen, and fly ash.

Material

- Body composition
Outer material: Spun-bond polypropylene
Filter media: Melt-blown filter media
Inner material: Spun-bond polypropylene
- Elastic Band composition
Polyester, Spandex
- Nose-piece composition
Aluminum

Dimensions

- Body Size
Length: 175 mm; Width: 89 mm
- Length of Nose-piece
118 mm
- Pleat Depth
15 mm

Textile Flammability

Class I flame spread

Sterility

Non-Sterile

Usage

Disposable, Single-Use only

Packaging

50 pcs/box, 10 box/cv
Total: 500 masks per case

Storage Conditions

Store under cool, clean, and dry conditions. Avoid excessive heat (over 40°C or 104°F)

Country of Origin

China

This technical data sheet is issued by A.R. Medicom Inc (Asia) Ltd